

nanoLC-1D™ plus system

The NanoLC that delivers the highest peptide detection sensitivity now features a dedicated rapid sample loading pump

The NanoLC-1D Plus system for proteomics includes all of the capabilities and benefits of the NanoLC-1D system—direct pumping for precise nanoscale gradients, and automated peak parking—with the added capability of a dedicated sample loading pump that can load large samples as rapidly as 30 $\mu\text{L}/\text{min}$. The NanoLC-1D Plus system is based on Eksigent's microfluidic flow control (MFC) technology, which uses continuous feedback to an electronically controlled pressure source to maintain precise nanoscale flow rates as low as 20 nL/min without flow splitting. Additionally, the system includes Eksigent's Dynamic Flow Control capability, which allows the flow rate to be changed rapidly for higher speed sample loading or for extended MS/MS analysis (peak parking) to identify co-eluting and low abundance peptides.

full integration and flexibility

The NanoLC-1D Plus system includes nanoscale binary gradient pumps, a dedicated sample loading pump, autosampler, and two injection valves. The 10-port injection valve accommodates 2 trap columns. The Eksigent Control Software provides complete control of the pumps, valves, and autosampler and provides flexible sample injection modes. The software is fully integrated in Thermo's Xcalibur® and ABI/MDS Sciex's Analyst® 1.4.1 software.

simple operation

The NanoLC-1D Plus system's intuitive control software simplifies setting up a rapid sample load followed by a nanoscale gradient. In addition to minimized mobile phase preparation and waste generation, direct pumping allows experimental conditions (e.g. column dimensions) to be changed easily without the need to reconfigure flow splitting hardware.

nanoLC-1D plus system

- High throughput—Rapid sample loading pump injects samples at up to 30 $\mu\text{L}/\text{min}$ onto trap or analytical columns, quickly loading large samples
- Precision—Direct pumping without flow splitting provides increased flow precision, resulting in improved protein identification
- Increased sensitivity—The NanoLC-1D Plus system incorporates Eksigent's unique Dynamic Flow Control capability, enabling peak parking, which improves the ability to identify low-abundance peptides
- Simple operation—the intuitive Eksigent Control Software includes full integration in Xcalibur and Analyst 1.4.1. Elimination of flow splitters eliminates time-consuming hardware reconfiguration





eksigent technologies
2021 las positas court, suite 161
livermore, california 94551
tel: 925 960 8869 fax: 925 960 8867
web: www.eksigent.com

new jersey office/laboratory
11 deer park, suite 204
monmouth junction, NJ 08852
tel: 732 274 9191

nanoLC-1D plus system specifications

configuration

Binary gradient pump, dedicated sample loading pump, 6-port low volume injection valve standard, 10-port optional.

gradient flow rate range

20 nL/min to 1000 nL/min (Other ranges available with flow module swap)

flow rate precision

Less than 0.5% RSD at 500 nL/min

sample loading pump

1 μ L/min to 30 μ L/min

pump type

Microfluidic flow control (MFC) direct pumping system with independent flow rate feedback for each mobile phase. Tunable PID flow rate control.

gradient formation

High-pressure gradient mixing, no flow splitting

delay volume

65 nL from mixer to external injection valve connected with 15 cm of 25 μ m id capillary.

delay time

<25 sec at 500 nL/min using 5 cm column and low port-to-port volume injection valve.

mobile phase compatibility

All mobile phases that are compatible with 316 stainless steel, PEEK, and silica.

injections

Standard and timed injections. Sample injection volume 200 nL to 50 μ L

columns

System optimized for use with 50 to 150 μ m id columns

peak parking

Peak parking flow rate user settable to as low as 5 nL/min. Triggered by user or through mass spectrometer.

peak parking response

Response time from 500 nL/min to 100 nL/min: <4 seconds at spray tip, using 150 μ m id, 6 cm long, 5 μ m C18 column

autosampler

One 96- or 384-well plate capacity, Peltier type temperature maintained sample compartment, 5 °C to 40 °C. Dual needle design for seal piercing. Integral 6-port injection valve. Numerous injection options.

system control

Computer with graphical user interface for control of all system parameters. Software creates data files for export and analysis. Drivers included for Xcalibur and Analyst 1.4.1

system size

14" (35 cm) wide x 8" (20 cm) high x 15" (38 cm) deep. Autosampler adds 18" (45 cm) in height.

utility requirements

100-240 VAC electrical input. Dry compressed air or nitrogen, 80-100 psi (6-7 bar)

GLP features

Tracking of instrument runtime, column usage, total injections, solvent usage, and error codes

diagnostics

System includes automated diagnostic routines for controller tuning, flow rate calibration, flow rate precision determination, and detection of leaks and blocks