

# Highly reproducible data from a microfluidic separation platform for application in proteomics

J. Bryce Young; Nicole Hebert; David Wyrick; Erika Lin; Remco van Soest  
 Eksigent Technologies 5875 Arnold Road, Dublin, CA 94568



## Overview

- Demonstration of the cHiPLC nanoLC platform as a high performance chromatographic tool for proteomics with the advantages of microfluidic design elements.
- Chromatographic performance on par with conventional packed capillary/spray tip combinations.
- High run-to-run reproducibility, and low inter-column variation.

## Introduction

- Variability in performance between nanoLC columns makes it difficult to standardize proteomic methods.
- Part of this variability is the highly skill dependent nature of making fluidic connections for nanoliter/min separations.
- Traditional nanoLC columns vary greatly, in packing bed length and quality.
- Microfluidic devices have inherent appeal as a chromatographic medium as the microchannels are defined lithographically and are highly reproducible.
- Microfluidic devices require a chip-to-world interface which does not sacrifice the performance gains from the microdevices themselves.
- In this poster, we will introduce a new microfluidic platform for nanoliter separations on microchips and demonstrate that robust, high quality separations enhances the reproducibility of proteomic experiments.

## Microfluidic Columns



The columns used in this work are microfabricated on fused silica using standard lithographic techniques.

- Circular channels in cross-section
- Packed at high pressure
- 3 μm particle packing

- Chip Column dimensions: 75 μm i.d. x 150 mm
- Packing material: Eksigent ChromXP C18 3μm 120Å

## Methods

- Separations were carried out on an Eksigent NanoLC Ultra 2Dplus LC system and a cHiPLC-nanoflex microfluidic system.




- A: 0.1% formic acid in H<sub>2</sub>O, and B: 0.1% formic acid in ACN.
- Gradient test method: 5 – 40%B at 2%/minute.
- Nanospray source from New Objective, using FS360-20-10, non-coated, pulled tips.
- MS data was collected using a Thermo LTQ ion trap.
- Data treatment and interpretation was performed using Mascot Distiller and MASCOT


- Samples in this work include:
  - Sigma Aldrich Peptide Mixture at 50 pg/ul
  - BSA, reduced, alkylated (carbamidomethyl), and digested with Trypsin

## Method Workflow

- cHiPLC nanoflex platform is operated in direct injection mode
  - Main valve is toggled in the "inject" position
  - Column eluate is plumbed directly into a nanospray source



## Results



- Represented to the left are 3 columns, each demonstrating a gradient separation of standard peptide mixture.
- Peaks 1 – 4 have retention times which are represented by RT1 – RT4.
- Inter-run variation is summarized for each column
- The table below summarizes the inter-column retention time reproducibility.

## Inter-column reproducibility

Inter-column	RT1	RT2	RT3	RT4
%STD DEV	0.60	1.26	0.76	0.82


- Inter-column reproducibility is key to experiments that run over long periods of time (months to years). e.g. validation of biomarkers.

## Inter-Column Sequence Coverage

- Two microfluidic columns are compared.
- Below is a Venn diagram that describes the inter-column overlap in sequence coverage

## Inter-Run Sequence Coverage

- A high degree of inter-run reproducibility is important in establishing robust proteomic methods.
- Below is a Venn diagram that describes the inter-run overlap in sequence coverage.



## Conclusions

- Microfluidic columns offer high inter-column reproducibility
- Inter-run repeatability is a measure of method robustness
- High quality, reproducible separations enable standardized methods for proteomics